DOC022.98.90169

# HT 200 S

10/2021, Edition 6

Basic User Manual

Basis Bedienungsanleitung

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## Specifications

Specifications are subject to change without notice.

HT 200 S	
<b>Type</b>	Fully automated high-temperature thermostat: 1000 watt high-performance heating Fast cooling through two ventilator system
<b>Digestion temperatures</b>	40–150 °C or 170 °C
<b>Digestion vessel compartments</b>	Twelve digestion compartments for HACH round cuvettes and reaction tubes
<b>Digestion times</b>	5–240 minutes
<b>Heating-up time</b>	Maximum 8 minutes
<b>Cooling-down time</b>	Maximum 13 minutes at an ambient temperature of < 25 °C
<b>Operating modes</b>	Three standard programs Nine possible user-specific programs
<b>Display</b>	Two-line alphanumeric display
<b>Temperature accuracy</b>	Thermostat block ± 3 °C
<b>Power consumption</b>	1300 VA
<b>Supply voltage</b>	230 VAC ± 10% 50 Hz
<b>Cold instrument socket fuse</b>	T 6.3A L; 250 V
<b>Dimensions</b>	Width 300 mm/height 330 mm/depth 430 mm (width 11.81 in./height 12.99 in./depth 16.93 in.)
<b>Weight</b>	10 kg (22 lb)

<b>HT 200 S</b>	
<b>Protection class</b>	I
<b>Pollution degree</b>	2
<b>Ovvervoltage category</b>	II
<b>Operating temperature</b>	16 to 29 °C (61 to 84 °F)
<b>Humidity</b>	maximum 80% relative humidity, non-condensing
<b>Altitude</b>	2000 m (6562 ft)
<b>Environmental conditions</b>	Indoor use

## General Information

### **WARNING**

Removing, bridging or disabling safety devices, safety functions and monitoring devices is forbidden.

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

## Safety information

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements.

Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

## Use of hazard information

### **DANGER**

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

### **WARNING**

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

### **CAUTION**

Indicates a potentially hazardous situation that may result in minor or moderate injury.

### **NOTICE**

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

## Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed.

This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.

This symbol indicates that a risk of electrical shock and/or electrocution exists.

This symbol indicates that the marked item can be hot and should not be touched without care.

Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems. Return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.

## Operating environment

### ⚠ WARNING

The instrument may not be used in dangerous environments.

The manufacturer and its suppliers reject any express or indirect guarantee for use with high-risk activities.

### NOTICE

Any use other than use in accordance with requirements defined in the user manual leads to the loss of the warranty claims and can lead to personal injury and property damage, for which the manufacturer assumes no liability.

Follow in addition to the local valid guidelines the following safety notes.

#### Safety notes for the regulatory use of the device:

- Do not operate the device near highly flammable substances, such as fuels, easily combustible chemicals and explosives.
- Do not operate the device near combustible gases, vapors or dust.
- Do not subject the device to strong vibrations or impacts.
- Do not open the device.
- Guarantee is voided if the instrument is not used in accordance with the guidelines present in this document.

Observe the following points to allow the instrument to function normally and give a long operating life.

- Position the instrument securely on a flat surface taking care to remove any objects from under the device.
- The ambient temperature for transport, storage and operation must be between +16 and +29 °C.

### NOTICE

Protect the instrument from extreme temperatures from heaters, direct sunlight and other heat sources.

- The relative humidity should be less than 80 %; moisture should not condense on the instrument.
- Leave at least a 15 cm clearance above and at each side of the instrument to allow air to circulate and to prevent electrical parts from overheating.
- Do not use or store the device in extremely dusty, humid or wet places.

#### Chemical and biological safety

### ⚠ DANGER

Potential danger with contact with chemical/biological substances.

Working with chemical samples, standards and reagents can be dangerous. Make yourself familiar with the necessary safety procedures and the correct handling of the chemicals before the work and read and follow all relevant safety data sheets.

Normal operation of this instrument may involve the use of hazardous chemicals or biologically harmful samples.

- Before handling these substances observe all danger notes and safety information printed on the containers of the original solutions and in the safety data sheet.
- Dispose of all consumed solutions in accordance with the national regulations and laws.

- Select the type of protective equipment according to the concentration and quantity of the dangerous substance at the respective work place.

## Start Up

### **WARNING**

The instrument may not be used in dangerous environments.

The manufacturer and its suppliers reject any express or indirect guarantee for use with high-risk activities.

### **NOTICE**

Use only an earthed socket for the connection of this device to the power supply. If you are not sure if the socket is earthed, have this checked by a qualified electrician.

The power plug serves in addition to the power supply to isolate the device quickly from the mains where necessary.

This is recommended for long-term storage and can prevent potential dangers in the event of a fault.

Therefore make sure that socket to which the device is connected is easy to reach by each user at all times.

### **NOTICE**

Any use other than use in accordance with requirements defined in the user manual leads to the loss of the warranty claims and can lead to personal injury and property damage, for which the manufacturer assumes no liability.

## Power connections

### **WARNING**

Electrical and fire hazards.

Use only the supplied power cable.

Only qualified experts may perform the tasks described in this section of the manual, while adhering to all locally valid safety regulations.

1. The power cable is connected to the rear of the thermostat with an earthed mains socket (230 volt +5 %/-15 %/50 Hz).
2. Switch on the instrument by pressing the toggle switch on the back of the instrument.

**Note:** *Do not turn the device off and on in rapid succession. Always wait approximately 20 seconds before switching on again so as not to damage the electronics and mechanics of the instrument.*

# Product overview

## Product contents

Please check that the order is complete. If anything is missing or damaged, please contact the manufacturer immediately.

The following components are supplied as standard with the thermostat:

- HT 200 S high-temperature thermostat
- Removal aid for cuvette adapter
- Cuvette adapter 20/13 mm (12 x)
- Power cable
- User Manual

**Note:** If any of these items are missing or damaged, please contact the manufacturer or the responsible sales representative immediately.


The HT 200 S high-temperature thermostat with HSD (high-speed digestion) technology meets all the technical requirements and safety regulations for fast and efficient digestion.

The HT 200 S high-temperature thermostat has twelve digestion compartments for round cuvettes and reaction tubes.

### NOTICE

Samples in 13-mm (LCK) round cuvettes cannot be digested unless the accompanying adapter is used! The adapter is removed from the thermostat using the removal tool; see figure 1.

Figure 1 Removal tool and adapter

1 Removal tool

2 Adapter

### DANGER

For quality and safety reasons, only HACH original cuvettes and reaction tubes may be used when conducting the analyses.

The sample cuvettes and sample reaction tubes are prepared in accordance with the currently valid working procedure.

Select the required temperature program on the HT 200 S in accordance with the information in the currently valid working procedure.


Higher temperatures than those specified in the working procedure can cause the cuvettes or reaction tubes to shatter or become damaged in the thermostat block.

If this happens, contact Technical Customer Service.

**Figure 2 Front view**

**Figure 3 Bottom of instrument**

1 Air circulation slots	5 LED – red, lights up when the lid locking mechanism is activated
2 Air intake vent	6 Display
3 Ventilation tube	7 Keys
4 Thermostat block	8 Air intake vent on bottom of the instrument

1 Air filter	3 Air filter grid
2 Filter mat	

## Thermostat block

The thermostat block is protected by the protective lid — even in the event of a cuvette or reaction tube smashing.

Any vapors and chemicals that may be released are contained in the thermostat block. They cannot penetrate inside the thermostat, so no subsequent technical defects can occur as a consequence of such releases.

### NOTICE

Damaged cuvettes and reaction tubes should not be used again (danger of glass breakage!).

Do not disconnect the power plug!

Organic solvents must not be heated in the instrument!

Only original HACH cuvettes and reaction tubes may be used!

Risk of injury. The locking pin on the housing lid has sharp edges

Thanks to the electrical locking feature, the thermostat offers the highest possible level of safety for the user. This locking feature for the digestion chamber closes automatically once the digestion process starts and only unlocks the lid once the digestion process is complete and the appropriate cooling-down phase has elapsed (exception: standard COD program). It is therefore impossible to open the thermostat during the digestion process.

Figure 4 Front view with lid open

1 Locking pin

## Switch on the instrument

1. Connect the power cable to the mains outlet.
2. Switch on the instrument by pressing the toggle switch on the back of the instrument.
3. The instrument starts automatically. The main menu is displayed after the current instrument version.

## Language selection

The thermostat software supports several languages. Press **key 1** when you switch on the thermostat.

Press **key 4** to scroll to the required language.

Then confirm the selection by pressing **key 3** under the **language symbol**.

The language is changed immediately and the start menu appears on the display.

## Main menu

### Overview of menu items

#### Digestion times

After selecting the **100** or **HT** standard program, the **Time** sub menu appears automatically; this sub menu is used to select the digestion duration/time.

#### Standard program 100

The digestion temperature is 100 °C. The **Time** sub menu then appears; this sub menu is used to select the required digestion time.

The following digestion times can be selected for a digestion temperature of 100 °C: 30, 60 and 120 minutes.

#### Standard program HT

The digestion temperature is 150 °C or 170 °C. The **Time** sub menu then appears; this sub menu is used to select the required digestion time.

The following digestion times can be selected at a digestion temperature of 150 °C or 170 °C: 10, 15 and 30 minutes.

The temperature of 170 °C shown on the display corresponds to the temperature of the heating block.

**Note:** The total digestion time consists of the selected digestion time and the heating-up and cooling-down times. This time can vary depending on the digestion temperature, the ambient temperature and the number of round cuvettes and/or reaction tubes.

At this heating block temperature, different temperatures are reached inside the digestion vessels, depending on the boiling point of the substance to be digested:


- Aqueous solutions containing diluted acids and lyes reach **150 °C** (e.g. LatoN, LCW 902 crack set, phosphate cuvette tests etc.).
- Semi-concentrated and concentrated acids (all COD cuvette tests) reach **170 °C**.

#### Standard program COD

The digestion temperature is **148 °C** and the digestion time is **120 minutes** (heating-up time included). In this program, it is not possible to select the digestion time.

COD digestion can be started immediately at 148 °C and 120 minutes.

Figure 5 Main menu — key assignment

**Note:** Wait at least 20 seconds before restarting the instrument in order to avoid damaging its electronics.

### Overview of additional menu items

Press **key 4** to display additional menu items **PRG** (Programming), **VAR** (Varying) and **CPRG** (Delete).

## PRG

Program user-specific digestions by entering the digestion term, temperature, digestion time and position

## VAR

Vary user-specific digestions by changing the digestion term, temperature, digestion time or position (= adjustment of the user-specific programmed parameters in the main menu)

## CPRG

Delete user-specific digestions.

Press **key 4** to return to the main menu.

## Start digestion

### DANGER

Potential danger with contact with chemical/biological substances.

Working with chemical samples, standards and reagents can be dangerous. Make yourself familiar with the necessary safety procedures and the correct handling of the chemicals before the work and read and follow all relevant safety data sheets.

The sample cuvettes and sample reaction tubes are prepared in accordance with the currently valid working procedure.

On the instrument, select the temperature program specified in the currently valid working procedure.

Higher temperatures than those specified in the working procedure can cause the cuvettes or reaction tubes to shatter or become damaged in the thermostat block. If this happens, contact Technical Customer Service.

### CAUTION

#### Burn hazard

Do not put a finger in an opening in the heating block.

Close the protective cover before increasing the temperature of the heating block, and always keep the protective cover closed during operation.

Burn hazard. The sample cuvettes are hot. Wear thermal insulating material (e.g. gloves or finger cots). Do not remove sample cuvettes if their temperature is above 100 °C.

Chemical hazards. If a reaction tube/cuvette breaks, do not let the liquid touch skin. If necessary, use a flue to convey chemical vapors.

### WARNING

The final cleaning of the thermostat should only be carried out by Technical Customer Service.

### NOTICE

To prevent reaction tube/cuvette and instrument damage, keep the openings in the heating block dry. Dry the outside of the reaction tubes/cuvettes fully.

To prevent damage to the instrument:

#### If glass breaks or fluid leaks out during digestion:

- 1 Press the **key** under **Esc** to abort the digestion. Let the instrument cool.
- 2 Turn the instrument off and disconnect the power cord.

#### If glass breaks or fluid leaks out when loading the thermostat:

- 1 Turn the instrument off and disconnect the power cord.
- 2 Remove the liquid with a pipette. Do not let the liquid touch skin.
- 3 Dispose of the liquid in accordance with regulations.

Then clean the instrument. See [section Clean the device, page 17](#).

Do not cover the ventilation slits in the protective cover. The instrument and reaction tubes/cuvettes may become too hot and decrease the accuracy of the measurement.

- 1 Prepare the cuvettes as described in the working procedure.

2. Seal the cuvette.
3. Carefully dry the outside of the cuvettes.
4. Place the cuvettes in the thermostat block and close the lid.

Once the digestion times have been set, the thermostat is ready for operation. The selected digestion conditions are shown in the top line of the display.


### Start the digestion in standard programs 100, HT or COD, or in a user-specific program

Press **key 1** under **Start** to start heating up the thermostat.

The lid is locked automatically (red LED lights up) at the same time (when the lid is closed).

**Note:** If the safety lid is open or incorrectly locked, the error message **Close lid** appears on the display after pressing **key 1** under **Start**. If the lid is then locked correctly after the error message has been displayed, the digestion starts automatically. The display is continuously updated once the digestion has started.

**Figure 6 Digestion display**

<b>1</b> <b>Digestion temperature</b>	<b>4</b> Once the <b>heating-up phase</b> is complete, the remaining digestion time is continuously updated and displayed (except for on the COD program: the digestion time starts being continuously updated during the heating-up phase).
<b>2</b> <b>Selected digestion time</b>	<b>5</b> These symbols flash during the heating-up phase. Once the <b>target temperature</b> has been reached, the current block temperature is displayed.
<b>3</b> <b>Abort the digestion</b> An <b>Abort</b> confirmation prompt appears. Press the <b>key</b> under <b>Yes</b> to abort the digestion; press the <b>key</b> under <b>No</b> to continue the digestion. If the heating-up phase is already finished when the digestion is aborted, there may be a cooling-down phase lasting a few minutes before the <b>Digestion display closes!</b> message appears on the display. Press the <b>key</b> under <b>Esc</b> to return to the <b>Time</b> sub menu.	

**Note:** When the standard program **COD** is selected, the digestion time is fixed at 120 minutes (including heating-up time).

When user-specific programs are used, the digestion time is specified by the user.

## Digestion end

### Standard programs 100 and HT, and user-specific digestion operations

The cooling-down phase (automatic rapid cooling) begins immediately after the digestion time has elapsed and after the digestion is aborted by pressing **key 4** under **Esc**. As soon as a temperature of **approx. 80 °C** is reached, the lock is automatically released; the **Close flap** instruction on the display must always be observed. The end of the cooling-down phase or the abortion of the digestion is signaled by an acoustic sound and by the **Digestion over!** message on the display.

Press **key 4** under **ESC** to return to the main menu.

### Standard program COD

The cooling-down phase is stopped at a temperature of **approx. 100 °C** and the lid lock is released. This is indicated by an acoustic signal.

The **Please invert cuvette** message appears on the display.

#### NOTICE

Take care when inverting — the cuvettes are hot! Wear protective clothing, protective gloves and safety glasses/face protection that are suitable for the work being performed.

After the cuvettes have been inverted and placed back in the thermostats, and the lid has been closed, **key 3** under the **arrow symbol** can be pressed to continue cooling to a temperature of approx. 44 °C. In this instance, the lid is locked again until the temperature has cooled to 80 °C.

The end of the cooling-down phase or the abortion of the digestion is indicated by an acoustic signal and by the **Digestion over!** message on the display.

Press **key 4** under **Esc** to return to the main menu.

## User programs

There are nine program locations available that can be configured by the user.

### Programming (PRG) of user-specific temperature programs

Figure 7 Programming

1	<b>Temperature</b> menu	4	Actuate key to select <b>VAR</b> . User-specific digestion programs are varied using the <b>VAR</b> menu.
2	Press key to return to the <b>Main menu</b> .	5	Press key to select <b>PRG</b> . User-specific digestions can be programmed using the <b>PRG</b> menu.
3	Press key to select <b>CPRG</b> . User-specific digestion programs can be deleted using the <b>CPRG</b> menu.		

User-specific digestion programs are programmed using the **PRG** menu. The **digestion term**, **temperature**, **digestion time** and **position** must be defined or calculated by the user.

### Enter the digestion term

Press **key 1** to set the saved alphanumeric character set to the required number, letter or space.

Press **key 2** to move the character shown as **\_** on the display one place to the right or back to the starting position.

A maximum of **three** characters can be entered.

Press **key 3** under **OK** to confirm the entered digestion term.

### Enter the digestion temperature

The digestion temperature is determined by the boiling point of the substance to be digested.

Press **key 1** to adjust the digestion temperature in 5 °C increments over a temperature range of 40 °C to 150 °C or 170 °C.

Press **key 3** under **OK** to confirm the entered digestion temperature.

### Influence of the acid concentration on the temperature reached in the sample containers

Selected digestion temperature (°C)	Temperature (°C) in the reaction tubes/ cuvettes for digestions of:		Thermostat block temperature (°C)
	Aqueous solutions with diluted acids and lyes	Semi-concentrated to concentrated acids	
40	40	41	41
60	60	62	62
80	80	84	84
100	100	110	110
120	120	134	134
130	130	146	146
135	135	152	152
140	140	158	158
150	150	170	170

**Note:** The programmed digestion temperature for the digestion of aqueous solutions with diluted acids and lyes (boiling points of 100 to 110 °C) is the temperature that is reached in the reaction tubes or cuvettes during the digestion process. The temperature in the thermostat block may be up to 20 °C higher!

For the digestion of semi-concentrated to concentrated acids (boiling

points of 150 °C and above), the temperatures reached in the reaction tubes or cuvettes are higher than the programmed digestion temperature.

In this instance, the temperatures in the reaction tubes or cuvettes are equivalent to the temperatures in the thermostat block.

### Enter the digestion time

Press **key 1** to adjust the digestion time (time) in 5 minute increments over a range of 5 to 240 minutes.

Press **key 3** under **OK** to confirm the entered digestion time.

### Position (first,, 2. position etc.) of the user-specific programs in the main menu

Press **key 1** to specify the required position (1–9) on the display.

Press **key 3** under **OK** to confirm the entered position.

### Open the user-specific digestion programs

Press **key 4** under the **arrow symbol** in the main menu to show additional menu items:

Press the key under the relevant digestion term to select the corresponding digestion program.

**Note:** This menu only opens if user-specific digestion programs have been programmed.

### Variation (VAR) of user-specific digestion programs

User-specific digestion programs are varied using the **VAR** menu.

Press **key 2** under **VAR** in the main menu.

The **Digestion term**, **Temperature**, Digestion time and Position parameters can be varied by the user.

### Select the digestion program to be varied

Press **key 1** under the **arrow symbol** in the main menu to show additional menu items.

Press the keys under the relevant digestion term to select the digestion program to be varied.

## Variation of the digestion term

Press **key 1** to set the saved alphanumeric character set to the required number, letter or space.

Press **key 2** to move the character shown as \_ on the display one place to the right or back to the starting position.

A maximum of **three** characters can be entered.

Press **key 3** under **OK** to confirm the changed digestion term.

## Variation of the digestion temperature

The digestion temperature is determined by the boiling point of the substance to be digested.

Press **key 1** to adjust the digestion temperature in 5 °C increments over a temperature range of 40 °C to 150 °C or 170 °C.

Press **key 3** under **OK** to confirm the changed temperature.

**Note:** The standard programs 100, HT and COD cannot be varied.

## Delete (CPRG) user-specific digestion programs

User-specific digestion programs are deleted using the **CPRG** menu.

Press **key 3** under **CPRG** in the main menu.

### Select the digestion program to be deleted


Press **key 2** under the **arrow symbol** in the main menu to show additional menu items.

Press the keys under the relevant digestion term to select the digestion program to be deleted.

Press **key 3** under **OK** to confirm the selection.

Press **key 3** under **OK** to delete the selected digestion program. The message will appear on the display

**Figure 8 Display — delete:**

1	Delete: message	2	XYZ? - Digestion term
3	Esc	4	Key 4
5	Key 3	6	Key 2
7	Key 1	8	OK

Press **key 1** under **OK** to delete the digestion program.

Press **key 4** under **Esc** to abort the deletion operation.

The display then automatically returns to the main menu.

**Note:** The standard programs 100, HT and COD cannot be deleted.

## Power off the instrument

Press the toggle switch on the back of the instrument.

## Error messages and display messages

Error message	Definition	Resolution
Prog.Mem. full	The user-specific program memory is full	Delete a digestion program that is no longer required
Prog.Mem. Clear	The user-specific program memory is empty	Enter a user-specific digestion program
Name not recognised	No name has been entered	Designation
Name already used esc	A name has been assigned twice	Select a different name
Close the lid	The lid has not been closed	Close the lid
ERR: 1000 temp. too high	Permissible temperature exceeded	Contact the Service department
ERR: 1001 Catch	Lock faulty	Contact the Service department

## Maintenance and service

Change the air filter (every year or if it is very dirty):

- Turn the instrument off and disconnect the power cord.
- Remove the fitted square air filter from the bottom of the instrument and replace the dirty filter mat with a new one.

### NOTICE

Check the condition of your air filter (level of dirt) on a regular basis, as a very dirty filter mat lengthens the cooling-down phase of the thermostat.

## Clean the device

### DANGER

Potential danger with contact with chemical/biological substances. Working with chemical samples, standards and reagents can be dangerous. Make yourself familiar with the necessary safety procedures and the correct handling of the chemicals before the work and read and follow all relevant safety data sheets.

### WARNING

Fire hazard! Do not use flammable agents to clean the instrument. To prevent instrument damage, do not use cleaning agents such as turpentine, acetone or similar products to clean the instrument or the display.

### CAUTION

Burn hazard  
Do not clean the instrument when it is hot.

The outside of the instrument can be cleaned with a dry cloth or a mild household cleaning agent. Under no circumstances may solvents such as petroleum spirit, acetone or similar substances be used.

- Turn the instrument off and disconnect the power cord.
- When the instrument is cool, clean the surface of the instrument with a soft, moist cloth and a weak soap solution. Make sure that no water gets into the instrument.

## Change the fuse in the cold instrument socket

### DANGER

Some circuits in this device are protected by fuses against overvoltages. For constant protection against risk of fire, replace these fuses only with fuses of the same type and power rating.

Defective fuses generally indicate a problem in the device. If fuse errors repeatedly occur, contact the service department for instructions on the return procedure for repair purposes. Under no circumstances attempt to repair the device independently.

### WARNING

Electrical and fire hazards.

Use only the supplied power cable.


Only qualified experts may perform the tasks described in this section of the manual, while adhering to all locally valid safety regulations.

1. Turn the instrument off and disconnect the power cord.
2. Wait until the temperature of the thermostat block and the sample cuvettes has dropped.
3. Remove the fuse holder by pressing the catch in the direction of the arrow.
4. Remove the faulty fuse (T6, 3A) (see [Figure 9](#))
5. Insert the accompanying spare fuse from the chamber.
6. Re-insert the fuse holder (clicks audibly into place)

### NOTICE

If the replacement fuse also blows after it has been inserted and the instrument has been switched on, contact our Technical Customer Service.

Figure 9 Back view

1	Air intake vent	5	Spare fuse
2	Protective grid	6	Fuse holder
3	Cold instrument socket	7	Toggle switch (on/off)
4	Fuse		

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)